

Anti-Ccdc6 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q16204

Item Number: CM0023735

Introduction

Rabbit polyclonal antibody against human CCDC6, validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for oncology and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	CCDC6 Rabbit pAb
Remarks/Alias	H4; PTC; TPC; PTC1; TST1; D10S170; CCDC6
Species	Human
GeneID	8030 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-222 of human CCDC6 (NP_005427.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16204
Protein Weight	53kDa
Shipping	Ice bag