


Anti-Ccdc61 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6R9

Item Number: CM0018695



Introduction

CCDC61 Rabbit Polyclonal Antibody validated for Western blot and ELISA in human samples. Research tool for centrosomal protein studies, centriole cohesion, and microtubule anchoring.

[Learn More](#)

Parameter	Value
Product Name	CCDC61 Rabbit pAb
Remarks/Alias	VFL3; hVFL3; CCDC61
Species	Human
Gene ID (Human)	729440
Gene ID	729440
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 352-531 of human CCDC61 (NP_001073871.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y6R9
Protein Weight	57kDa
Shipping	Ice bag