

Anti-Ccdc88A (Girdin) Polyclonal Antibody For Wb, Elisa - Q3V6T2

Item Number: CM0022035



Introduction

High-quality anti-CCDC88A (Girdin) polyclonal antibody validated for Western blotting (WB) and ELISA, targeting human CCDC88A/Girdin (UniProt Q3V6T2). Ideal for cell migration, EGFR signaling, and cytoskeleton studies.

[Learn More](#)

Parameter	Value
Product Name	CCDC88A Rabbit pAb
Remarks/Alias	APE; GIV; GRDN; PEHO; HkRP1; PEHOL; GIRDIN; KIAA1212; CCDC88A
Species	Human
GeneID (Human)	55704
GeneID	55704
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1571-1870 of human CCDC88A (NP_001129069.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q3V6T2
Protein Weight	216kDa
Shipping	Ice bag