

Anti-Ccdc92 Rabbit Polyclonal Antibody For Wb, Elisa - Q53Hc0

Item Number: CM0030663

Introduction

Rabbit polyclonal antibody targeting human CCDC92, validated for Western blotting and ELISA. CCDC92 plays a role in innate immunity and inhibits ebolavirus replication. Suitable for research use.

[Learn More](#)

Parameter	Value
Product Name	CCDC92 Rabbit pAb
Remarks/Alias	CCDC92
Species	Human
GeneID (Human)	80212
GeneID	80212
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q53HC0
Protein Weight	37kDa
Shipping	Ice bag