

Anti-Cck Rabbit Polyclonal Antibody For Wb And Elisa - P06307

Item Number: CM0019942



Introduction

Rabbit polyclonal antibody against human CCK, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studies of gall bladder contraction, pancreatic enzyme release, and gastric acid secretion.

[Learn More](#)

Parameter	Value
Product Name	CCK Rabbit pAb
Remarks/Alias	CCK
Species	Human
GeneID (human)	885
GeneID	885
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-115 of human CCK (NP_000720.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P06307
Protein Weight	13kDa
Shipping	Ice bag