

Anti-Cckar (Cck1-R) Rabbit Polyclonal Antibody For Wb And Elisa - P32238

Item Number: CM0022558



Introduction

CCKAR rabbit polyclonal antibody for WB and ELISA. Detects human CCKAR, cross-reacts with mouse. Ideal for studying cholecystokinin receptor type A in feeding, gall bladder contraction, and signaling.

[Learn More](#)

Parameter	Value
Product Name	CCKAR Rabbit pAb
Remarks/Alias	CCK-A; CCK1R; CCKRA; CCK1-R; CCKAR
Species	Human
Gene ID	886 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 100-200 of human CCKAR (NP_000721.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P32238
Protein Weight	48kDa
Shipping	Ice bag