

Anti-Cckbr Rabbit Polyclonal Antibody For Wb, Elisa - P32239

Item Number: CM0025809



Introduction

Rabbit polyclonal antibody targeting human CCKBR (GASR/CCK-B/CCK2R). Recognizes Human, Mouse, Rat. Validated for WB and ELISA. Ideal for gastrin/CCK signaling, neurobiology, and cancer research. - P32239

[Learn More](#)

Parameter	Value
Product Name	CCKBR Rabbit pAb
Remarks/Alias	GASR; CCK-B; CCK2R; CCKBR
Species	Human
GeneID (Human)	887
GeneID	887
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-335 of human CCKBR (NP_795344.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32239
Protein Weight	48kDa
Shipping	Ice bag