


Anti-Ccl1 (I-309) Rabbit Polyclonal Antibody For Wb, Elisa - P22362

Item Number: CM0013319



Introduction

CamelBio supplies rabbit anti-human CCL1 (I-309) polyclonal antibody for WB and ELISA. This chemokine attracts monocytes via CCR8, serving as a key tool for inflammation and immune response research.

[Learn More](#)

Parameter	Value
Product Name	CCL1 Rabbit pAb
Remarks/Alias	P500, S1Se, TCA3, I-309, SCYA1, CCL1
Species	Human
GeneID (Human)	6346
GeneID	6346
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-96 of human CCL1 (NP_002972.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P22362
Protein Weight	11kDa
Shipping	Ice bag