

Anti-Ccl17 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q92583

Item Number: CM0010335



Introduction

Rabbit monoclonal antibody against human CCL17 (TARC, SCYA17). Validated for WB, IHC-P, ELISA with cross-reactivity to rat. Ideal for studying Th2-mediated inflammation, allergy, and chemokine signaling.

[Learn More](#)

Parameter	Specification
Product Name	CCL17 Rabbit mAb
Remarks/Alias	TARC; ABCD-2; SCYA17; A-152E5.3
Species	Human
Gene ID (Human)	6361
Gene ID	6361
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-94 of human CCL17 (NP_002978.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q92583
Protein Weight	11kDa
Shipping	Ice bag