


Anti-Ccl19 Polyclonal Antibody For Wb, If-P, Elisa - Q99731

Item Number: CM0020077



Introduction

CamelBio offers rabbit polyclonal anti-CCL19 antibody validated for WB, IF-P, ELISA. Recognizes human, mouse, rat CCL19, a key chemokine for T/B-cell trafficking and inflammation research - Q99731.

[Learn More](#)

Parameter	Value
Product Name	CCL19 Rabbit pAb
Remarks/Alias	ELC; CKb11; MIP3B; MIP-3b; SCYA19; CCL19
Species	Human
GeneID (Human)	6363
GeneID	6363
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-95 of human CCL19 (NP_006265.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99731
Protein Weight	11kDa
Shipping	Ice bag