

Anti-Ccl22 Rabbit Polyclonal Antibody For Ihc-P, Elisa - 000626

Item Number: CM0024371

Introduction

CamelBio offers CCL22 Rabbit Polyclonal Antibody (000626) validated for IHC-P and ELISA. Detects human CCL22 (MDC) with cross-reactivity to rat. Ideal for immune trafficking and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CCL22 Rabbit pAb
Remarks/Alias	MDC; ABCD-1; SCYA22; STCP-1; DC/B-CK; A-152E5.1; CCL22
Species	Human
GeneID (Human)	6367
GeneID	6367
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-93 of human CCL22 (NP_002981.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	O00626
Protein Weight	11kDa
Shipping	Ice bag