


Anti-Ccl28 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrj3

Item Number: CM0026974



Introduction

Rabbit polyclonal antibody targeting human CCL28 chemokine, validated for WB and ELISA. Detects human, mouse, and rat CCL28; involved in T-cell and eosinophil chemotaxis. For research use in mucosal immunity studies.

[Learn More](#)

Parameter	Value
Product Name	CCL28 Rabbit pAb
Remarks/Alias	MEC; CCK1; SCYA28; CCL28
Species	Human
GeneID	56477
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-127 of human CCL28 (NP_683513.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRJ3
Protein Weight	14kDa
Shipping	Ice bag