

Anti-Ccl2/Mcp-1 Polyclonal Antibody For Wb And Elisa - P13500

Item Number: CM0024651



Introduction

Rabbit polyclonal antibody against human CCL2/MCP-1, validated for WB and ELISA. Target is a monocyte chemoattractant involved in inflammation and atherosclerosis. Ideal for chemokine research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CCL2/MCP-1 Rabbit pAb
Remarks/Alias	HC11; MCAF; MCP1; MCP-1; SCYA2; GDCF-2; SMC-CF; HSMCR30; CCL2/MCP-1
Species	Human
Gene ID (Human)	6347
Gene ID	6347
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-99 of human CCL2/MCP-1 (NP_002973.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13500
Protein Weight	11kDa
Shipping	Ice bag