

Anti-Ccl2/Mcp-1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P13500

Item Number: CM0029133



Introduction

Rabbit polyclonal antibody against human CCL2/MCP-1, validated for WB, IHC-P, and ELISA. Detects endogenous levels of total CCL2/MCP-1 with cross-reactivity to mouse and rat. Ideal for inflammation and atherosclerosis research.

[Learn More](#)

Parameter	Value
Product Name	CCL2/MCP-1 Rabbit pAb
Remarks/Alias	HC11, MCAF, MCP1, MCP-1, SCYA2, GDCF-2, SMC-CF, HSMCR30, CCL2/MCP-1
Species	Human
GeneID (Human)	6347
GeneID	6347
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-99 of human CCL2/MCP-1 (NP_002973.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13500
Protein Weight	11kDa
Shipping	Ice bag