

Anti-Ccl3 (Mip-1-Alpha) Rabbit Polyclonal Antibody For Wb, Elisa - P10147

Item Number: CM0029669



Introduction

High-quality rabbit polyclonal antibody against human CCL3 (MIP-1-alpha, SCYA3) for WB and ELISA. Suitable for inflammation, chemokine, and HIV-suppressive factor research.

[Learn More](#)

Parameter	Value
Product Name	CCL3 Rabbit pAb
Remarks/Alias	MIP1A; SCYA3; G0S19-1; LD78ALPHA; MIP-1-alpha; CCL3
Species	Human
GeneID (Human)	6348
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-92 of human CCL3 (NP_002974.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10147
Protein Weight	10kDa
Shipping	Ice bag