

Anti-Ccl3 Rabbit Polyclonal Antibody For Wb, Elisa - P10147

Item Number: CM0029283



Introduction

High-quality rabbit polyclonal antibody targeting human CCL3 (MIP-1-alpha), validated for Western blotting and ELISA applications with cross-reactivity to mouse and rat. Ideal for inflammatory and HIV research.

[Learn More](#)

Parameter	Value
Product Name	CCL3 Rabbit pAb
Remarks/Alias	MIP1A; SCYA3; G0S19-1; LD78ALPHA; MIP-1-alpha; CCL3
Species	Human
Gene ID (Human)	6348
Gene ID	6348
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 24-92 of human CCL3 (NP_002974.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10147
Protein Weight	10kDa
Shipping	Ice bag