

Anti-Ccl3L1 Rabbit Polyclonal Antibody For Wb, Elisa - P16619

Item Number: CM0029383

Introduction

Rabbit polyclonal antibody against human C-C motif chemokine 3-like 1 (CCL3L1, SWISS P16619). Validated for WB and ELISA, suitable for studying chemotaxis, HIV inhibition, and inflammatory diseases.

[Learn More](#)

Parameter	Value
Product Name	CCL3L1 Rabbit pAb
Remarks/Alias	LD78; 464.2; MIP1AP; SCYA3L; G0S19-2; SCYA3L1; D17S1718; LD78BETA; LD78-beta(1-70); CCL3L1
Species	Human
Gene ID (Human)	6349
Gene ID	6349
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16619
Protein Weight	10kDa
Shipping	Ice bag