


Anti-Ccl4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P13236

Item Number: CM0022145



Introduction

This anti-CCL4 rabbit polyclonal antibody detects human CCL4 (MIP-1 β), a chemokine involved in inflammation and HIV suppression. Suitable for WB, IF/ICC, and ELISA. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	CCL4 Rabbit pAb
Remarks/Alias	ACT2; G-26; HC21; LAG1; LAG-1; MIP1B; SCYA2; SCYA4; MIP1B1; AT744.1; MIP-1-beta; CCL4
Species	Human
GeneID (human)	6351
GeneID	6351
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-92 of human CCL4 (NP_002975.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P13236
Protein Weight	10 kDa
Shipping	Ice bag