

Anti-Ccl4/Mip-1 Beta Monoclonal Antibody For Wb, If/Icc, Elisa - P14097

Item Number: CM0005212

Introduction

High-quality rabbit monoclonal antibody against mouse CCL4/MIP-1 beta (P14097). Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for inflammation research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	CCL4/MIP-1 beta Rabbit mAb
Remarks/Alias	Act-2; Mip1b; Scya4; MIP-1B; AT744.1
Species	Mouse
GenelD (Human)	6351
GenelD	20303
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14097
Protein Weight	10 kDa
Shipping	Ice bag