

Anti-Ccl5 Polyclonal Antibody For Wb, Elisa - P50231

Item Number: CM0020509

Introduction

Rabbit polyclonal antibody against rat CCL5 (C-C motif chemokine 5, RANTES). Validated for WB and ELISA. Ideal for chemokine signaling and inflammatory response research.

[Learn More](#)

Parameter	Value
Product Name	CCL5 Rabbit pAb
Remarks/Alias	C-C motif chemokine 5
Species	Rat
Gene ID (Human)	6352
Gene ID	81780
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-52 of rat CCL5 (AAA96499.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P50231
Protein Weight	10kDa
Shipping	Ice bag