

Anti-Ccl5/Rantes Polyclonal Antibody For Wb, If-P, Elisa - P13501

Item Number: CM0016605



Introduction

Rabbit polyclonal antibody against human CCL5/RANTES, validated for WB, IF-P, ELISA with cross-reactivity to mouse and rat. Suitable for chemokine research and HIV inhibition studies.

[Learn More](#)

Parameter	Value
Product Name	CCL5/RANTES Rabbit pAb
Remarks/Alias	SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta; CCL5/RANTES
Species	Human
GeneID (Human)	6352
GeneID	6352
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-91 of human CCL5/RANTES (NP_002976.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13501
Protein Weight	10kDa
Shipping	Ice bag