

Anti-Ccl8 Rabbit Polyclonal Antibody For Ihc-P And Elisa - P80075

Item Number: CM0029436



Introduction

Rabbit polyclonal antibody targeting human CCL8 (MCP-2), suitable for IHC-P and ELISA. Detects human and mouse CCL8, an 11 kDa chemokine involved in monocyte chemotaxis and inflammation.

[Learn More](#)

Parameter	Value
Product Name	CCL8 Rabbit pAb
Remarks/Alias	HC14; MCP2; MCP-2; SCYA8; SCYA10; CCL8
Species	Human
GeneID (Human)	6355
GeneID	6355
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-99 of human CCL8 (NP_005614.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P80075
Protein Weight	11kDa
Shipping	Ice bag