


Anti-Ccl8/Mcp-2 Rabbit Polyclonal Antibody For Western Blot - P80075

Item Number: CM0021107



Introduction

Rabbit polyclonal antibody targeting human CCL8/MCP-2, suitable for Western blot. Detects the 11kDa chemotactic cytokine that attracts monocytes, lymphocytes, and eosinophils. For research use.

[Learn More](#)

Parameter	Value
Product Name	CCL8/MCP-2 Rabbit pAb
Remarks/Alias	CCL8; HC14; MCP-2; MCP2; SCYA10; SCYA8; C-C motif chemokine 8; CCL8/MCP-2
Species	Human
GeneID (Human)	6355
GeneID	6355
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-99 of human CCL8/MCP-2 (NP_005614.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P80075
Protein Weight	11kDa
Shipping	Ice bag