

Anti-Ccm2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bsqs5

Item Number: CM0028613



Introduction

Rabbit polyclonal antibody targeting CCM2 (malcavernin), a key regulator of cardiovascular development. Suitable for WB, IHC-P, ELISA, with reported cross-reactivity to human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	CCM2 Rabbit pAb
Remarks/Alias	OSM; C7orf22; PP10187; CCM2
Species	Human
Gene ID	83605
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CCM2 (NP_113631.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BSQ5
Protein Weight	49kDa
Shipping	Ice bag