

Anti-Ccnb1Ip1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npc3

Item Number: CM0019169



Introduction

Rabbit polyclonal antibody against human CCNB1IP1, validated for Western blot and ELISA. Cross-reacts with mouse. Ideal for investigating meiotic crossover, cell cycle regulation, and ubiquitin ligase activity.

[Learn More](#)

Parameter	Value
Product Name	CCNB1IP1 Rabbit pAb
Remarks/Alias	HEI10; C14orf18; CCNB1IP1
Species	Human
Gene ID (Human)	57820
Gene ID	57820
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-50 of human CCNB1IP1 (NP_878269.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NPC3
Protein Weight	32kDa
Shipping	Ice bag