


Anti-Ccr2 Rabbit Polyclonal Antibody For Wb And Elisa - P41597

Item Number: CM0027430



Introduction

CCR2 rabbit polyclonal antibody validated for WB and ELISA, reactive with human and mouse. Ideal for chemokine receptor signaling, monocyte migration, and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	CCR2 Rabbit pAb
Remarks/Alias	CKR2; CCR-2; CCR2A; CCR2B; CD192; CKR2A; CKR2B; CMKBR2; MCP-1-R; CC-CKR-2; CCR2
Species	Human
GeneID	729230
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human CCR2 (NP_001116513.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P41597
Protein Weight	42kDa
Shipping	Ice bag