

Anti-Ccr3 Rabbit Polyclonal Antibody (Ckr3/Cd193) For Wb, Ihc-P, Elisa - P51677

Item Number: CM0020237



Introduction

High-quality rabbit anti-CCR3 polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat CCR3. Recombinant immunogen; detects ~41kDa band. Ideal for chemokine receptor studies in immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	CCR3 Rabbit pAb
Remarks/Alias	CKR3; CD193; CKR 3; CMKBR3; C C CKR3; CC-CKR-3; CCR3
Species	Human
GeneID (Human)	1232
GeneID	1232
Immunogen	Recombinant fusion protein containing a sequence within amino acids 146-226 of human CCR3 (NP_001828.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51677
Protein Weight	41kDa
Shipping	Ice bag