

Anti-Ccr5 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P51681

Item Number: CM0023018



Introduction

Anti-CCR5 rabbit polyclonal antibody validated for WB, IF/ICC, IF-P, and ELISA. Detects human, mouse, and rat CCR5. Ideal for studying chemokine signaling, HIV entry, and inflammatory diseases.

[Learn More](#)

Parameter	Value
Product Name	CCR5 Rabbit pAb
Remarks/Alias	CKR5; CCR-5; CD195; CKR-5; CCCKR5; CMKBR5; IDDM22; CC-CKR-5; CCR5
Species	Human
GeneID (Human)	1234
GeneID	1234
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 246-345 of human CCR5 (NP_000570.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51681
Protein Weight	41kDa
Shipping	Ice bag