

Anti-Ccr6/Cd196 Rabbit Monoclonal Antibody (Fc, Elisa) - O54689

Item Number: CM0005624

Introduction

Recombinant rabbit monoclonal antibody against mouse CCR6/CD196, for flow cytometry and ELISA. Ideal for research on chemokine receptor signaling, T-cell trafficking, and mucosal immunity in inflammatory and autoimmune disease models.

[Learn More](#)

Parameter	Value
Product Name	CCR6/CD196 Rabbit mAb
Remarks/Alias	CCR-6; KY411; Cmkbr6; CC-CKR-6
Species	Mouse
Gene ID (Human)	1235
Gene ID	12458
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-38 of mouse CCR6/CD196 (NP_033965.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Mouse
SWISS	O54689
Protein Weight	42kDa
Shipping	Ice bag