


Anti-Ccr7 Rabbit Monoclonal Antibody For Wb, Elisa - P32248

Item Number: CM0002336



Introduction

High-quality Anti-CCR7 (CD197) rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat CCR7, a chemokine receptor critical for lymphocyte migration. Ideal for immunology and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	CCR7 Rabbit mAb
Remarks/Alias	BLR2; EBI1; CCR-7; CD197; CDw197; CMKBR7; CC-CKR-7; CCR7
Species	Human
GeneID (Human)	1236
GeneID	1236
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CCR7 (P32248).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32248
Protein Weight	43kDa
Shipping	Ice bag