


Anti-Ccr7 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P32248

Item Number: CM0015503



Introduction

CCR7 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA. It recognizes human, mouse, and rat CCR7 (C-C chemokine receptor type 7), with a predicted MW of 43kDa. Ideal for immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CCR7 Rabbit pAb
Remarks/Alias	BLR2; EB11; CCR-7; CD197; CDw197; CMKBR7; CC-CKR-7; CCR7
Species	Human
Gene ID (Human)	1236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 330-372 of human CCR7 (NP_001829.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32248
Protein Weight	43kDa
Shipping	Ice bag