

Anti-Ccr8 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P51685

Item Number: CM0005249



Introduction

High-quality rabbit monoclonal antibody targeting human CCR8, validated for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for chemokine receptor research in immunology and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	CCR8 Rabbit mAb
Remarks/Alias	CY6; TER1; CCR-8; CKRL1; CDw198; CMKBR8; GPRCY6; CMKBRL2; CC-CKR-8; CCR8
Species	Human
GeneID (human)	1237
GeneID	1237
Immunogen	A synthetic peptide corresponding to amino acids 250-355 of human CCR8 (P51685).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51685
Protein Weight	41kDa
Shipping	Ice bag