


Anti-Ccr9 Rabbit Monoclonal Antibody For Wb, Elisa - P51686

Item Number: CM0010037



Introduction

Anti-CCR9 Rabbit Monoclonal Antibody validated for WB and ELISA. Reacts with human, mouse, rat CCR9 (42kDa). Suitable for chemokine receptor and HIV-1 coreceptor studies.

[Learn More](#)

Parameter	Value
Product Name	CCR9 Rabbit mAb
Remarks/Alias	GPR28; CDw199; GPR-9-6; CC-CKR-9; CCR9
Species	Human
Gene ID (Human)	10803
Gene ID	10803
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CCR9 (NP_112477.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51686
Protein Weight	42kDa
Shipping	Ice bag