


Anti-Ccr9 Rabbit Polyclonal Antibody For Wb, Elisa - P51686

Item Number: CM0013353



Introduction

This rabbit polyclonal antibody targets C-C chemokine receptor type 9 (CCR9), a key mediator of immune cell trafficking. Validated for WB and ELISA, it recognizes human, mouse, and rat CCR9, making it suitable for immunology and HIV studies.

[Learn More](#)

Parameter	Value
Product Name	CCR9 Rabbit pAb
Remarks/Alias	GPR28; CDw199; GPR-9-6; CC-CKR-9; CCR9
Species	Human
Gene ID (Human)	10803
Gene ID	10803
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 300-369 of human CCR9 (NP_112477.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51686
Protein Weight	42kDa
Shipping	Ice bag