

Anti-Ccs Rabbit Polyclonal Antibody For Wb And Elisa - O14618

Item Number: CM0025079



Introduction

Rabbit polyclonal antibody against human CCS (copper chaperone for superoxide dismutase), suitable for Western blot and ELISA. Recombinant immunogen corresponding to amino acids 1-274. For research use.

[Learn More](#)

Parameter	Value
Product Name	CCS Rabbit pAb
Remarks/Alias	CCS
Species	Human
GeneID (Human)	9973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human CCS (NP_005116.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14618
Protein Weight	29kDa
Shipping	Ice bag