


Anti-Cct3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P49368

Item Number: CM0029026



Introduction

CCT3 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Detects human, mouse, rat CCT3, a TRiC chaperonin subunit involved in actin, tubulin folding and ciliogenesis. Ideal for chaperone research.

[Learn More](#)

Parameter	Value
Product Name	CCT3 Rabbit pAb
Remarks/Alias	CCTG; PIG48; TRIC5; CCT-gamma; TCP-1-gamma; CCT3
Species	Human
Gene ID	7203 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CCT3 (NP_005989.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49368
Protein Weight	61kDa
Shipping	Ice bag