


Anti-Cct4 Polyclonal Antibody For Wb, If/Icc, Elisa - P50991

Item Number: CM0027398



Introduction

CCT4 Rabbit pAb is a polyclonal antibody validated for WB, IF/ICC, and ELISA, targeting human CCT4 with cross-reactivity to mouse and rat. Suitable for studying chaperonin function in protein folding, telomere maintenance, and ciliogenesis.

[Learn More](#)

Parameter	Value
Product Name	CCT4 Rabbit pAb
Remarks/Alias	SRB; Cctd; CCT-DELTA; CCT4
Species	Human
GeneID (Human)	10575
GeneID	10575
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 130-390 of human CCT4 (NP_006421.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50991
Protein Weight	58kDa
Shipping	Ice bag