


Anti-Cct5 Polyclonal Antibody For Wb, Elisa - P48643

Item Number: CM0026245



Introduction

CCT5 rabbit polyclonal antibody validated for Western blot and ELISA applications. Recognizes human, mouse, and rat CCT5, a key component of the TRiC chaperonin complex assisting actin and tubulin folding.

[Learn More](#)

Parameter	Value
Product Name	CCT5 Rabbit pAb
Remarks/Alias	CCTE; HEL-S-69; PNAS-102; CCT-epsilon; TCP-1-epsilon; CCT5
Species	Human
GeneID (Human)	22948
GeneID	22948
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human CCT5 (NP_036205.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48643
Protein Weight	60kDa
Shipping	Ice bag