


Anti-Cct6A Polyclonal Antibody For Wb, Elisa - P40227

Item Number: CM0026883



Introduction

CCT6A Rabbit pAb (P40227) is a polyclonal antibody validated for WB and ELISA. It detects human, mouse, and rat CCT6A, a TRiC chaperonin subunit that assists actin and tubulin folding.

[Learn More](#)

Parameter	Value
Product Name	CCT6A Rabbit pAb
Remarks/Alias	CCT6; Cctz; HTR3; TCPZ; TCP20; MoDP-2; TTCP20; CCT-zeta; CCT-zeta-1; TCP-1-zeta; CCT6A
Species	Human
GeneID (human)	908
GeneID	908
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-250 of human CCT6A (NP_001753.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40227
Protein Weight	58kDa
Shipping	Ice bag