


Anti-Cct6A Rabbit Polyclonal Antibody For Wb And Elisa - P40227

Item Number: CM0021110



Introduction

CCT6A Rabbit Polyclonal Antibody validated for Western blot and ELISA. Recognizes human and rat CCT6A, a chaperonin subunit involved in actin and tubulin folding. Suitable for proteostasis research.

[Learn More](#)

Parameter	Specification
Product Name	CCT6A Rabbit pAb
Remarks/Alias	CCT6; Cctz; HTR3; TCPZ; TCP20; MoDP-2; TTCP20; CCT-zeta; CCT-zeta-1; TCP-1-zeta
Species	Human
Gene ID (Human)	908
Gene ID	908
Immunogen	Recombinant protein corresponding to a sequence within amino acids 80-250 of human CCT6A (NP_001753.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P40227
Protein Weight	58kDa
Shipping	Ice bag