


Anti-Cct7 Polyclonal Antibody For Wb, Elisa - Q99832

Item Number: CM0016261



Introduction

Rabbit polyclonal antibody targeting human CCT7 (CCT-eta, TCP1 subunit eta), validated for Western blot and ELISA. Detects endogenous levels in human, mouse, and rat. Ideal for chaperonin and protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	CCT7 Rabbit pAb
Remarks/Alias	CCTH; CCTETA; NIP7-1; TCP1ETA; CCT7
Species	Human
GeneID (human)	10574
GeneID	10574
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-500 of human CCT7 (NP_006420.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99832
Protein Weight	59kDa
Shipping	Ice bag