


Anti-Cct8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P50990

Item Number: CM0026868



Introduction

Rabbit anti-CCT8 polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat CCT8. Suitable for protein folding, telomere biology, and cilia research.

[Learn More](#)

Parameter	Value
Product Name	CCT8 Rabbit pAb
Remarks/Alias	Cctq; PRED71; D21S246; C21orf112; CCT8
Species	Human
GeneID (Human)	10694
GeneID	10694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 248-497 of human CCT8 (NP_001269837.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50990
Protein Weight	60kDa
Shipping	Ice bag