

Anti-Cd10 Rabbit Monoclonal Antibody () For Flow Cytometry - P08473

Item Number: CM0008358



Introduction

-conjugated rabbit monoclonal antibody targeting human CD10 (Neprilysin/CALLA). Validated for flow cytometry, enabling detection of hematopoietic malignancies and normal lymphoid subsets. Ideal for immunophenotyping and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD10 mAb
Remarks/Alias	NEP; SFE; CD10; CALLA; CMT2T; SCA43
Species	Human
Gene ID (Human)	4311
Gene ID	4311
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of human CD10 (NP_000893.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08473
Protein Weight	86kDa
Shipping	Ice bag