

Anti-Cd10 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - Q61391

Item Number: CM0005457



Introduction

Rabbit monoclonal antibody specific to mouse CD10 (Neprilysin), validated in WB, FC, ELISA. Supports diagnostic development and research on peptide metabolism and skin remodeling.

[Learn More](#)

Parameter	Value
Product Name	CD10 Rabbit mAb
Remarks/Alias	Mme
Species	Mouse
GeneID (Human)	4311
GeneID	17380
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-750 of mouse CD10 (NP_032630.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Mouse
SWISS	Q61391
Protein Weight	86kDa
Shipping	Ice bag