

Anti-Cd101 Polyclonal Antibody For Wb, Elisa - Q93033

Item Number: CM0022951

Introduction

High-quality rabbit polyclonal antibody against human CD101 (Immunoglobulin superfamily member 2), validated for Western blot and ELISA. Useful for investigating T-cell proliferation inhibition and leukemic tumor cell research.

[Learn More](#)

Parameter	Value
Product Name	CD101 Rabbit pAb
Remarks/Alias	V7; IGSF2; EWI-101; CD101
Species	Human
GeneID (human)	9398
GeneID	9398
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101 (NP_004249.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag