

Anti-Cd101 Rabbit Monoclonal Antibody For Wb And Elisa - Q93033

Item Number: CM0004315

Introduction

Rabbit monoclonal antibody against human CD101 (V7, IGSF2, EWI-101) for Western blot and ELISA. Inhibits T-cell proliferation; leukemic marker. Immunogen: recombinant protein. Swiss-Prot Q93033.

[Learn More](#)

Parameter	Value
Product Name	CD101 Rabbit mAb
Remarks/Alias	V7; IGSF2; EWI-101; CD101
Species	Human
Gene ID	9398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101. (NP_004249.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag