

Anti-Cd101/Igsf2 Monoclonal Antibody For If/Icc And Flow Cytometry - Q93033

Item Number: CM0008803



Introduction

Rabbit monoclonal antibody against human CD101 (IGSF2, V7, EWI-101). Validated for IF/ICC and flow cytometry. Detects a ~115 kDa membrane glycoprotein; suitable for T-cell inhibition and leukemic marker research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD101/IGSF2 mAb
Remarks/Alias	V7; IGSF2; EWI-101
Species	Human
Gene ID (Human)	9398
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101/IGSF2 (NP_004249.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	Q93033
Protein Weight	115kDa
Shipping	Ice bag