

## Anti-Cd101/Igsf2 Rabbit Monoclonal Antibody For Elisa, Wb, If/Icc, Fc - Q93033

Item Number: CM0007418

### Introduction

Rabbit monoclonal antibody against human CD101/IGSF2 (V7, EWI-101). Validated for ELISA, WB, IF/ICC, FC. Cross-reacts with mouse. Suitable for T-cell proliferation and immunology research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD101/IGSF2 Rabbit mAb                                                                                                  |
| Remarks/Alias    | V7; IGSF2; EWI-101                                                                                                      |
| Species          | Human                                                                                                                   |
| GenelD           | 9398 (Human)                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 21-410 of human CD101/IGSF2 (NP_004249.2) |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | ELISA, WB, IF/ICC, FC                                                                                                   |
| Cross Reactivity | Human, Mouse                                                                                                            |
| SWISS            | Q93033                                                                                                                  |
| Protein Weight   | 115kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                 |