

Anti-Cd102 Rabbit Polyclonal Antibody For Wb, Elisa - P35330

Item Number: CM0019126



Introduction

Rabbit polyclonal antibody against mouse CD102 (ICAM-2). Validated for Western blot and ELISA. Reacts with mouse and rat. Target is a cell adhesion molecule involved in leukocyte migration and immune responses.

[Learn More](#)

Parameter	Value
Product Name	CD102 Rabbit pAb
Remarks/Alias	CD102; Ly-60; Icam-2
Species	Mouse
GeneID (Human)	3384
GeneID (Mouse)	15896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-222 of mouse CD102 (NP_034624.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35330
Protein Weight	31kDa
Shipping	Ice bag