


Anti-Cd102 Rabbit Polymab® For Wb, Ihc-P, Elisa - P13598

Item Number: CM0007246



Introduction

CD102 Rabbit PolymAb® suitable for WB, IHC-P, ELISA. Recognizes human, mouse, rat ICAM2, a 30 kDa cell adhesion molecule important for immune response.

[Learn More](#)

Parameter	Specification
Product Name	CD102 Rabbit PolymAb®
Remarks/Alias	CD102
Species	Human/Mouse
GeneID (Human)	3384
GeneID	3384/15896
Immunogen	Recombinant protein of human CD102
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13598
Protein Weight	30KDa
Shipping	Ice bag